

NexTemp® Plus Single-Use Clinical Thermometer with Rectal Sheaths

2117

PRODUCT SPECIFICATIONS

Storage:

At or Below 86°F (30°C)

Supplier FDA Registration ID:

2246308

Product Information:

Lead Time: 0-14 days

Minimum Order Requirement: 1 Case

Packaging: 50 ct. box / 40 boxes per case

Shipping Information:

	Box	Case
Weight (lbs.)	0.15	7.00
Dimensions (in.)	0.75 x 4.5 x 4.75	10.25 x 10.25 x 8.75

MEDICAL INDICATORS NexTemp® Plus
 REF 2117-20
 2000 50 40 1
 LOT 99999 2999-12 CE 0086 Latex Free Celsius 38°C 98°F
 (01) 601 694 4840 1052 (17) 2999-12 (10) 99999
 Made in U.S.A. by Medical Indicators Inc, Hamilton, NJ 08619 15-2117-20-04

LOT 99999 2999-12
 (01) 30169448401051 (17) 2999-12 (10) 99999
 15-2117-01 99999

ABOUT THE PRODUCT

NexTemp® Plus Clinical Thermometers measures oral or axillary temperature in Celsius Scale.

Each dot contains a different, precise mixture that changes color at a specific temperature. Activated dots change from metallic green to black with the last black dot indicating the patient's temperature.

Temperature Range: 35.5-40.4°C

Accuracy: +/- 0.1°C (37.0°C through 39.0°C)
 +/- 0.2°C (elsewhere)

- Latex-Free
- Non-Sterile



50 Thermometers with
 50 Rectal Sheaths Per Box

MEDICAL INDICATORS

Medical Indicators, Inc.
 16 Thomas J. Rhodes Industrial Drive
 Hamilton, NJ 08619
 Office: (609) 737-1600
 Fax: (609) 587-8635
 info@medicalindicators.com